**VANGUARD STEERABLE DIAGNOSTIC CATHETER ID BWAB**  
**VANGUARD STEERABLE DIAGNOSTIC CATHETER BI BWAC**  
**VANGUARD STEERABLE DIAGNOSTIC CATHETER BI ID BWAD**GEBRAUCHSANWEISUNG  
INSTRUCTIONS FOR USE  
MODE D'EMPLOI  
GEBRUIKSAANWIJZING  
MANUAL DE INSTRUCCIONES  
INSTRUÇÕES DE UTILIZAÇÃO  
BRUKSANVISNING  
BRUKSANVISNING  
UPUTE ZA UPOTREBUDeutsch  
English  
Français  
Nederlands  
Español  
Português  
Norsk  
Svensk  
Hrvatski Vanguard AG  
Landsberger Straße 266  
12623 Berlin, Deutschland  
www.vanguard.de  
service@vanguard.de

©Vanguard AG 100001910\_REV\_B\_10-05-2019

**Verpackungsbeschriftung und –symbole / Packaging and labelling symbols**  
**/ Étiquetage et symboles sur l'emballage / Verpakkingstekst en –symbolen**  
**/ Rotulación y símbolos del envase / Rótulo e símbolos da embalagem /**  
**Påskrift og symboler på forpakningen / Text och symboler på förpackningen /**  
**Natpisi i simboli na ambalaži**

Inhalt / Content / Contenu / Inhoud / Contenido / Conteúdo / Innhold / Innehåll / Sadržaj:

Ein / one / un / een / un / um / én / en / jedan:

Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD

Medizinprodukt / Medical device / Dispositif médical / Medisch hulpmiddel / Producto médico / Dispositivo médico / Medisinsk utstyr / Medicinteknisk produkt / Medicinski proizvod

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## EN

### GENERAL

Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at [service@vanguard.de](mailto:service@vanguard.de) or contact our field staff before using the product on patients.

### General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

✓ Requirements

→ Instructions

• Lists

⚠ Safety instructions

The original version of these instructions for use is in German.

### Device Information

This product may only be operated if its safe use is guaranteed. Observe the handling and safety instructions!

Vanguard assumes no liability for damages that arise from improper operation of the Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD and/or disregard of these instructions for use.

Vanguard AG and the manufacturer of the original products, Biosense Webster, Inc. (33 Technology Drive, Irvine, California 92618, USA) are not affiliated companies; remanufacturing of the product is carried out solely by Vanguard AG, which places it on the market without the cooperation of Biosense Webster, Inc.

„Biosense Webster“ and „CARTO“ are trademarks of Biosense Webster, Inc. or one of its affiliated companies.

### PRODUCT DESCRIPTION

The Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD is equipped with a torsion-resistant shaft and a deflectable tip with several platinum electrodes. This catheter is intended for temporary intracardiac recording of potentials and temporary stimulation of the heart during electrophysiological examinations. Ablations may not be performed with the catheter.

The deflection of the tip (curvature) is controlled via a movable thumb knob (model ID BWAB) or a rocker lever (models bi BWAC / bi ID BWAD) at the handle respectively. The thumb knob is pushed forward and the rocker lever is turned respectively to bend the tip. The tip is straight when the thumb knob is pulled back and the rocker lever is in neutral position. The shape of the curvature depends on the length of the deflectable tip and the position of the thumb knob or rocker lever respectively. There are different catheters with different curvatures available (see "Technical Data"). The torsion resistant shaft allows the plane of the curvature to be rotated. In this way, the catheter can be positioned at the desired site. An extension cable with a suitable connector is used to connect to a suitable standard recording device.

The Vanguard Steerable Diagnostic Catheter ID BWAB / bi ID BWAD have a programmable memory device (EEPROM) that stores information for unique catheter identification. CARTO® EP Navigation Systems with Auto ID Technology are able to read this information and automatically recognize the catheter. These catheters have to be used with appropriate interface cables. Refer to appropriate Carto® EP Navigation System user manual for details on system settings related to the catheter.

### INTENDED PURPOSE

#### Indication

The Vanguard Steerable Diagnostic Catheter BWAB / bi BWAC / bi ID BWAD is intended for temporary intracardiac electrophysiological recording and stimulation (mapping) during the evaluation of cardiac arrhythmias. The catheter is designed for use in the coronary sinus.

#### Guidelines for use

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purpose about the safety and effective use of the Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The medical device may only be used by a physician who is appropriately trained in electrophysiology and familiar with its techniques. All electrophysiological examinations must be performed in a fully equipped electrophysiology laboratory/OR.

All Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD have the same compatibility with accessories (e.g. cables) as the corresponding products (see technical data) of the original Biosense Webster, Inc.

#### Contraindications

The application of the Vanguard Steerable Diagnostic Catheter ID BWAB | bi BWAC | bi ID BWAD is contraindicated for the following:

• Patients with intracardiac parietal thrombus

• Patients who have undergone a ventriculostomy or atriotomy in the preceding four weeks

• Patients with an aortic valve prosthesis,

• Transeptal access for patients with thrombus or myxoma in the left ventricle or with an intracardiac patch or baffle,

• Patients who have shown sensitivity to foreign objects in the past or who have an allergy

• Patients with severe anomalies in the heart and/or blood vessels (transposition) and obstruction

An active systemic infection is considered to be a relative contraindication; therefore, an examination should only be carried out after a rigorous risk/benefit assessment by the treating physician.

#### Side effects

When using the catheter, breakage of the catheter, dislocation of the catheter, contact problems between the catheter and pacemaker and an increase in the stimulus threshold can cause the following undesired side effects:

• pain and discomfort in the chest

• increase in phosphokinase levels

• subcutaneous hematoma

• stimulation of skeletal muscles and nerves

• lesion of the brachial plexus

• lesion of the thoracic duct

• lesion of the intracardiac conduction system

• damage to heart valves and/or blood vessels (e.g. venous supplies)

- vascular bleeding
- thrombosis
- vasovagal reactions
- local and systemic infection
- pneumothorax due to subclavian puncture
- haematothorax
- pericarditis
- cardiac tamponade
- arteriovenous fistula
- pseudo aneurysms
- arrhythmias (e.g. ventricular tachyarrhythmia)
- cardiac thromboembolism, air embolisms
- perforation of the atrium or ventricle (myocardium/endocardium) with or without loss of simulation or ventricular fibrillation
- coronary artery spasm, coronary artery thrombosis, coronary artery transection
- transient ischemic attack (TIA) and apoplectic stroke
- myocardial infarction
- stroke
- total heart block
- death

Previously unknown side effects may occur at any time. The distributing company must be notified of any side effects that are not specified.

In connection with the known side effects, please also refer to the relevant specialist literature.

## TECHNICAL DATA

### Vanguard Steerable Diagnostic Catheter ID BWAB

<b>Usable Length:</b>	115 cm
<b>Shaft diameter:</b>	6 F = 2 mm
<b>Number of electrodes:</b>	10
<b>Electrode spacing:</b>	2-8-2 mm

Item No.	Original ref. no. <sup>1</sup>	Curve type	Dimensions tip electrode / ring electrodes	Connection <sup>2</sup>	Additional Function
35087	D135303	D-Typ	2 mm / 1,3 mm	12 Pin; Redel	Auto ID
35050	D135304	F-Typ	2 mm / 1,3 mm	12 Pin; Redel	Auto ID

### Vanguard Steerable Diagnostic Catheter bi BWAC/ bi ID BWAD

<b>Usable Length:</b>	115 cm
<b>Shaft diameter:</b>	7 F = 2,3 mm
<b>Number of electrodes:</b>	10
<b>Electrode spacing:</b>	2-8-2 mm

Item No.	Original ref. no. <sup>1</sup>	Curve type	Dimensions tip electrode / ring electrodes	Connection <sup>2</sup>	Additional Function
35051	36D-35R	D-F -Typ	2 mm / 0,8 mm	10 Pin; Redel	-
35052	36G-35QCT	D-F -Typ	2 mm / 0,8 mm	12 Pin; Redel	Auto ID
35088	36G-5JMCT	F-J -Typ	2 mm / 0,8 mm	12 Pin; Redel	Auto ID

Other data → see label.

<sup>1</sup> Corresponding REF. No. of the original manufacturer

<sup>2</sup> "Redel" is a protected brand of INTERLEMO HOLDINGS S.A., Ecublens, CH or the brand of a manufacturer associated with this company.

This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC and all applicable changes.

## SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, consequence and counter-measure. A safety instruction indicates a potential risk to the health or life of persons.

⚠ Excessive bending or kinking of the catheter!

**Liquid penetration and impairment with the function of the catheter, damage to the wires and/or injury of the patient**

→ Prevent kinking of the catheter shaft or excessive bending.

→ Do not bend the distal tip manually.

⚠ Proximal plug connection or handle comes into contact with liquid!

**Interference with the function of the catheter**

→ Prevent the connector from coming into contact with liquids.

⚠ Plug connector damaged!

**Interference with electrical safety and function**

→ Carry out a careful visual inspection and functional inspection of the plug.

⚠ Electrical contacts of the catheter or connected cables are electrically connected to the heart!

**Ventricular fibrillation due to uncontrolled electrical energy**

→ Do not touch the plug contacts of the catheter or of connected cables with your fingers or with objects.

→ Only connect the catheter's plug contacts to medical devices using connection cables which meet the requirements and are intended for this purpose.

→ Do not allow the catheter's plug or the connection cable's plug to come into contact with liquids. This can impair the electrical safety and functionality.

⚠ Electromagnetic interference due to an external pacemaker!

**Damage or interference of other devices, ventricular fibrillation and burns**

→ Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.

⚠ Transfer of electrostatic discharges to the catheter system!

**Ventricular fibrillation due to uncontrolled electrical energy**

→ Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.

⚠ Long fluoroscopy times with X-rays!

**Somatic and genetic damage due to high exposure to ionizing radiation**

→ Particularly for pregnant women, children, young patients and patients who have to undergo several examinations, reconsider the risk/benefit ratio as a consequence of the potential effects of radiation.

⚠ Sterilized package is damaged or opened!

**Infection of the patient due to unsterile products**

→ Check the packaging before opening it and, if necessary, do not use the catheter.

⚠ Damage due to improper handling and/or unsuitable storage conditions!

**Damage to the central circulatory system, infection, septicemia**

→ Take note of the intended use and the relevant identification on the packaging.

→ Do not expose the catheter to any organic solvents such as alcohol.

⚠ Improper catheter use or the thumb knob not pushed back!

**Injuries, perforations and tamponades in the heart and vascular system**

→ Ensure that only trained personnel uses the catheter.

→ Always pull the thumb knob backwards to insert or remove the catheter in order to put the catheter tip into a straight position.

→ If there is resistance, do not use force to push the catheter forwards or backwards.

⚠ Reprocessing, including resterilization and reuse of the catheter by the user or operator!

**Injury of blood vessels and intracardiac injuries, infection, septicemia**

→ Ensure that the catheter is not reprocessed, reesterilized or reused.

## HANDLING INSTRUCTIONS

The catheter must be handled with great care in order to avoid injuries, perforations and tamponades in the heart and vascular system.

### Storage

The product must be stored in a cool, dry and light-protected area. Storage temperatures from 5 °C to 25 °C.

→ The information on storage conditions is also provided on the label of the product packaging.

### Safety measures before use

• Before using the catheter or external pacemaker, safety must read and understand the respective instructions for use.

• When connecting external devices and using them for their specific purpose, consult the corresponding instructions for use. Only use authorized medical devices under consideration of the respective instructions for use.

• All static electricity must be kept away from the catheter system. You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.

• A prerequisite for making an intracardiac potential recording and/or stimulation is the appropriate training of medical staff and the availability of a fully equipped electrophysiology laboratory. The procedure itself may only be carried out by a trained, experienced physician.

### Removal from packaging

→ Check that the packaging is not damaged. If the packaging is damaged the product has to be deposited and another new product for the procedure is has to be used.

→ Remove the catheter under aseptic conditions and transfer it to a sterile working environment.

→ Ensure that the packaging material is disposed of safely.

→ Check the catheter for any obvious defects such as kinks in the shaft, loose electrodes or damaged plugs and, if necessary, swap it for another catheter. The catheter must not be used if there is any obvious damage.

→ Carry out a careful visual inspection and functional inspection of the plug and replace the catheter if necessary.

### Application

✓ In order to avoid injuring the ventricular system and the heart, the intervention must always be carried out under X-ray control.

✓ An external pacemaker and a defibrillator must be available and ready for use throughout the entire procedure.

✓ The principles of asepsis must be observed.

→ Apply the standard method for vessel punctures, introduction of a guide wire and the application of an introducer sheath in compliance with the instructions for use.

→ Connect the catheter to a suitable recording device or Carto® EP Navigation System with an appropriate connection cable.

→ Before inserting and removing the catheter, ensure that the deflectable tip is straightened.

→ Guide and monitor the catheter movements within the patient with X-ray imaging and electrocardiograms.

→ Take care that the electrodes have a stable contact with the mapping site.

**NOTE:** In order to make positioning of the tip of Vanguard Steerable Diagnostic Catheter ID BWAB easier, it can be bent using the thumb knob. To do this, push the thumb knob forward (distal direction). If you push the thumb knob backwards, the tip will straighten itself again. The catheter is straight when the thumb knob is pulled back fully and can then be removed.

**NOTE:** With the models Vanguard Steerable Diagnostic Catheter bi BWAC/ bi ID BWAD with bidirectional steering the curvature can be bent to two opposed curves by rotating the rocker lever. The tip is straightened when the rocker lever is in neutral position (vertical to handle). The friction control knob allows the physician to adjust the resistance for moving the deflectable tip. Rotating the friction control knob in clockwise direction increases the amount of force; rotation in counter-clockwise direction decreases it.

**NOTE:** If resistance occurs, the catheter must not be pulled back with force. In this case, an X-ray must be used to check the position of the catheter.

### Disposal

Any used products must be considered as contaminated hospital waste. The used catheters and the protective covers must be disposed of in accordance with the applicable national laws and guideline.

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