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Read these instructions carefully before using the product described. If you have questions concerning the instructions for use or handling of the product, please contact our Customer Service Centre at service@vanguard.de or contact our field staff before using the product on patients.

General information

These instructions contain certain notations which are intended to help you understand the functions and meaning of the text more quickly:

✓	Requirements
→	Instructions
•	Lists
⚠	Safety instructions

All serious incidents related to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The original version of these instructions for use is in German.

Device information

This product may only be operated if its safe use is ensured by the user by following the instructions. Follow the notes and safety precautions!

Vanguard assumes no liability for damages that arise from improper operation of the Vanguard Steerable Diagnostic Catheter SJA | SJF and/or disregard of these instructions for use.

Vanguard AG and the manufacturer of the original product, *Abbott Medical (formerly Irvine Biomedical, Inc. a St. Jude Medical Company, 2375 Morse Avenue, Irvine, California 92614 USA)*, are not affiliated companies; the device is remanufactured solely by Vanguard AG, who markets the device without collaboration by the *Abbott Medical (formerly St. Jude Medical, Inc.)*.

DEVICE DESCRIPTION

The Vanguard Steerable Diagnostic Catheter SJA | SJF is a catheter equipped with a torsion-resistant shaft and a steerable tip with multiple platinum electrodes. This catheter is designed for temporary intracardiac recording of potentials and temporary stimulation of the heart. The catheter may not be used for ablation.

The catheter has the following clinical benefit: It is used as part of the technical setup for electrophysiological studies, together with recording / stimulation devices, further diagnostic and ablation catheters, electro anatomical mapping systems, and imaging devices. The device helps identifying cardiac structures that can be targets for catheter-based ablation. It supports acute procedural success, which is defined as successful restoration of sinus rhythm in atrial arrhythmias (AF / SVT) and elimination of inducible ventricular tachycardia in ventricular arrhythmias. The long-term clinical benefit offered by acutely successful procedures is reduced probability of arrhythmia recurrence, thus avoidance/reduction of cardiovascular risks/symptoms and improvement of quality of life. Based on literature data, an acute procedural success rate of at least 90% for atrial fibrillation, at least 75% for various types of supraventricular tachycardia, and at least 50% for ventricular tachycardias and premature ventricular contractions represents current state of the art for cardiac electrophysiology procedures.

A sliding piston with thumb knob on the handle controls the deflection of the tip (curvature). To deflect the tip, the thumb knob is pushed forward or the rocker lever is turned. If the thumb knob is in the rear position the tip is straight. The shape of the curve depends on the length of the steerable tip and the position of the thumb knob. Catheters with different curvatures are available (see "Technical data" section). The tip shape designations correspond to those of the original manufacturer *Abbott Medical (formerly St. Jude Medical, Inc.)*. The curvature plane can be rotated by turning the

torsion-resistant shaft. This allows the catheter to be positioned as desired. For connection to a suitable standard registration device, a connection cable with a suitable plug connection must be used.

Vanguard Steerable Diagnostic Catheter SJA

Usable length: 110 cm

REF	Corresponding REF of the original manufacturer	Shaft diameter	Curvature	Electrodes			Connector (corresponds to St.Jude Cable)
				Quantity	Spacing [mm]	Width [mm]	
35951	IBI-81532	4F	Large	10	2-5-2	1	14 Pin (1910)
36096	IBI-81534	4F	Large	10	5	1	14 Pin (1910)
35952	IBI-81540	4F	Medium	4	2-5-2	1	9 Pin (1904)
35949	IBI-81530	4F	Medium	10	2	1	14 Pin (1910)
35950	IBI-81531	4F	Medium	10	2-5-2	1	14 Pin (1910)
35948	IBI-81474	5F	Large	4	2-5-2	1	9 Pin (1904)
35938	IBI-81174	5F	Large	10	2-5-2	1	14 Pin (1910)
35947	IBI-81472	5F	Medium	4	2-5-2	1	9 Pin (1904)
35937	IBI-81172	5F	Medium	10	2-5-2	1	14 Pin (1910)
35944	IBI-81404	6F	Large	4	2-5-2	2	9 Pin (1904)
35945	IBI-81405	6F	Large	4	5	2	9 Pin (1904)
35955	IBI-81807	6F	Large	8	2	1	14 Pin (1910)
35956	IBI-81809	6F	Large	8	2-5-2	1	14 Pin (1910)
35933	IBI-81104	6F	Large	10	2-5-2	1	14 Pin (1910)
35935	IBI-81107	6F	Large	10	5	1	14 Pin (1910)
35936	IBI-81108	6F	Large	10	2	1	14 Pin (1910)
35942	IBI-81402	6F	Medium	4	2-5-2	2	9 Pin (1904)
35943	IBI-81403	6F	Medium	4	5	2	9 Pin (1904)
35953	IBI-81801	6F	Medium	8	2	1	14 Pin (1910)
35954	IBI-81802	6F	Medium	8	2-5-2	1	14 Pin (1910)
35932	IBI-81102	6F	Medium	10	2-5-2	1	14 Pin (1910)
35946	IBI-81418	6F	X-Large	4	2-5-2	2	9 Pin (1904)
35934	IBI-81105	6F	X-Large	10	2-5-2	1	14 Pin (1910)

REF	Corresponding REF of the original manufacturer	Shaft diameter	Curvature	Electrodes			Connector (corresponds to St.Jude Cable)
				Quantity	Spacing [mm]	Width [mm]	
35940	IBI-81207	7F	Super Large	20	5	1	26 Pin (1924)
35941	IBI-81209	7F	Super Large	20	2-5-2	1	26 Pin (1924)
35939	IBI-81202	7F	XX-Large	20	2-10-2	1	26 Pin (1924)

Other data → see label.

Vanguard Steerable Diagnostic Catheter SJF

Usable length: 110 cm

REF	Corresponding REF of the original manufacturer	Shaft diameter	Curvature	Electrodes			Connector (corresponds to St.Jude Cable)
				Quantity	Spacing [mm]	Width [mm]	
35928	IBI-81734	5F	Large	10	2-5-2	1	14 Pin (1910)
35929	IBI-81735	5F	Large	10	5	1	14 Pin (1910)
35927	IBI-81730	5F	Medium/ Large	10	2-8-2	1	14 Pin (1910)
35930	IBI-81945	6F	Large	10	2-5-2	1	14 Pin (1910)
35931	IBI-81947	6F	Medium/ Large	10	5	1	14 Pin (1910)

Other data → see label.

INTENDED PURPOSE

Operating principle

The Vanguard Steerable Diagnostic Catheter SJA | SJF is intended for temporary intracardiac electrogram recording, electrical stimulation, and mapping of the heart within the scope of electrophysiological examinations.

Indication

The Vanguard Steerable Diagnostic Catheter SJA | SJF is intended to be used together with compatible external hardware for electrophysiological studies of the heart ("EP studies"). Such EP studies are conducted to confirm, differentiate, or treat cardiac arrhythmias. Depending on indication, the medical device may be used within the atria, the ventricles, or the coronary sinus.

Contraindications

The application of the Vanguard Steerable Diagnostic Catheter SJA | SJF is contraindicated for the following:

- intracardiac parietal thrombus in patients who have undergone a ventriculostomy or atriotomy in the preceding four weeks,
- mechanical prosthetic heart valve through which the catheter must be inserted,
- echocardiographically confirmed and visible thrombus or myxoma,

- interatrial closure in persistent foramen ovale if transseptal access is required,
- aortic valve prosthesis with retrograde transaortic approach,
- sensitivity to foreign objects or allergies,
- severe anomalies of the heart and/or blood vessels (transposition) and obstruction,
- acute systemic infection,
- if catheter manipulation is unsafe (e.g. intracardiac mural thrombus),
- an embolic vena cava protective filter device and/or known femoral thrombus requiring femoral catheter insertion,
- for patients in whom a Vanguard Steerable Diagnostic Catheter SJA | SJF is inserted and magnetic resonance imaging (MRI) is required.

Target population and user group

The Vanguard Steerable Diagnostic Catheter SJA | SJF will be used by trained cardiologists specialized in electrophysiology or fellows in training for adult patients with suspected or confirmed cardiac arrhythmias and indication for an electrophysiological study as determined by the treating physician.

Use environment

The Vanguard Steerable Diagnostic Catheter SJA | SJF must only be used within a fully-equipped operational EP lab.

RESIDUAL RISKS

Side effects/potential adverse events

Besides risks that may be caused by inappropriate handling / lack of caution when operating the catheter electrophysiological studies generally are associated with complication risks independent of the specific equipment used. Type and frequency of complications is partly dependent on the intended procedure. The following typical risks have been identified as expected side effects by indication from literature review:

Atrial Fibrillation

Frequent (>1%)

- Access site complications (e.g. hematoma, arterio-venous fistula, pseudoaneurysm)
- Pericardial complications (hemorrhage, effusion, tamponade)

Occasional (0.1% - 1%)

- Cardiac valve damage
- Cerebrovascular embolism (stroke, transitory ischemic attack)
- Death
- Embolism
- Esophageal injury / atrio-esophageal fistula
- Infection (e.g. pericarditis, myocarditis, pneumonia, urinary tract infections)
- Phrenic nerve injury / palsy
- Pneumothorax / hemothorax
- Pulmonary vein stenosis

Supraventricular tachycardia

Occasional (0.1% - 1%)

- Access site complications (e.g. hematoma, arterio-venous fistula, pseudoaneurysm)
- Atrio-ventricular block
- Cerebrovascular embolism (stroke, transitory ischemic attack)
- Death
- Infection (e.g. pericarditis, myocarditis, pneumonia, urinary tract infections)
- Pericardial complications (hemorrhage, effusion, tamponade)

Rare (< 0.1%)

- Cardiac valve damage
- Embolism
- Myocardial infarction
- Phrenic nerve injury / palsy
- Pneumothorax / hemothorax

Ventricular tachycardia

Frequent (>1%)

- Access site complications (e.g. hematoma, arterio-venous fistula, pseudoneurysm)
- Atrio-ventricular block
- Cerebrovascular embolism (stroke, transitory ischemic attack)
- Coronary artery spasm / stenosis
- Death
- Embolism
- Heart failure / decompensation
- Incessant VT/VF
- Myocardial infarction
- Pericardial complications (hemorrhage, effusion, tamponade)

Occasional (0.1% - 1%)

- Bundle branch block
- Cardiac valve damage
- Cardiogenic shock
- Infection (e.g. pericarditis, myocarditis, pneumonia, urinary tract infections)
- Pneumothorax / hemothorax
- Pleural cavity puncture

Rare (<0.1%)

- Phrenic nerve injury / palsy

Previously unknown side effects may occur at any time. The user and patient must report all serious incidents associated with the medical device to the manufacturer and the appropriate authority. It is the responsibility of the user to inform the patient of the listed contraindications and adverse reactions associated with the use of a Vanguard Steerable Diagnostic Catheter SJA | SJF.

SAFETY INSTRUCTIONS

The safety instructions consist of three aspects: Cause, Consequence, Countermeasure. A safety instruction indicates a potential risk to the health or life of persons.

Electrical contacts of the catheter or connected cables are electrically connected to the heart!

Ventricular fibrillation due to uncontrolled electrical energy

- Do not touch the plug contacts of the catheter or of connected cables with your fingers or with objects.
- Only connect the catheter's plug contacts to medical devices using connection cables which meet requirements and are intended for this purpose.
- Do not allow the catheter's plug or the connection cable's plug to come into contact with liquids.

 Proximal plug connection or handle or cable comes into contact with liquid!

Impairment of the function of the catheter

- Prevent the plug, the cable and the handle from coming into contact with liquids.
- Do not immerse the plug and/or handle in solutions.
- Ensure that the cable and catheter connection remains dry throughout the procedure.

 Do not touch the plug contacts with your fingers or thin objects!

Impairment of the function of the catheter

- Ensure that the plug contacts do not come into contact with thin objects. Do not touch the plug contacts.

 The medical personnel has not been adequately trained!

Fault during intracardiac potential recordings and/or stimulation

- It is a pre-requirement for intracardiac potential recording and/or stimulation that specialist medical personnel are trained sufficiently and that the procedure is performed in a fully equipped EP laboratory. The procedure itself may only be carried out by a trained cardiologist specialised in electrophysiology.

 Insufficient caution when handling the catheter!

Injuries of the heart, perforation or tamponades

- Take great care when handling the catheter.

 Ignition of flammable gases or other substances!

Fire or explosion

- Ensure that flammable and easily inflammable materials are removed from the working area.

 Catheter is damaged!

Patient and/or user injury

- Subject the catheter to a careful visual and functional check.
- Damaged catheters may not be used.
- If necessary, replace damaged products with new replacement products.

 Electromagnetic interference due to an external pacemaker!

Damage or impairment of other devices/ventricular fibrillation and burns

- Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.
- All recording and stimulation devices must be electrically shielded.
- Equipment used with the Vanguard Steerable Diagnostic Catheter SJA | SJF must meet the following requirements: Type CF, defibrillation resistant, electrical safety design according to EN 60601-1 as well as all official requirements for the respective intended use.

 Transfer of electrostatic discharges to the catheter system!

Ventricular fibrillation due to uncontrolled electrical energy

- Only use the catheter in a properly equipped and operated medical electrophysiology laboratory.
- All recording and stimulation devices must be electrically insulated.
- Equipment used with the Vanguard Steerable Diagnostic Catheter SJA | SJF must meet the following requirements: Type CF, defibrillation resistant, electrical safety design according to EN 60601-1 as well as all official requirements for the respective intended use.

 Misplacement of the catheter into the coronary system during transaortic access!

Increased risk of bleeding, shock or heart attack

- Ensure that there is good imaging during transaortic access.

 Exposure to long fluoroscopy times with X-rays for patients and laboratory staff!

Somatic and genetic damage due to high exposure to ionising radiation

- The potential radiation exposure must be considered for the procedure and appropriate measures taken to minimise this exposure.
- Particularly for pregnant women, children, young patients and patients who have to undergo several examinations, reconsider the risk/benefit ratio as a consequence of the potential effects of radiation.

 Sterilised package is damaged or opened or the "use by" date is exceeded!

Infection of the patient due to non-sterile products

- Check the packaging before opening it. The catheter may not be used if its sterile packaging is open or damaged.
- The catheter may no longer be used after its "use by" date has expired.

 Damage due to continued use contrary to the intended use and due to unsuitable storage conditions!

Damage to the central circulatory system, infection, septicaemia

- Take note of the intended use and the relevant storage requirements on the packaging.

 Improper use or incorrect operation of the curvature control!

Injuries, perforations and tamponades in the heart and vascular system

- Ensure that only trained personnel uses the catheter.
- Vanguard Steerable Diagnostic Catheter SJA | SJF: When inserting or removing the catheter, always push the thumb knob backwards to straighten the catheter tip.
- If there is resistance, do not use force to push the catheter forwards or backwards.
- The catheter must not be subjected to any mechanical stress caused by manual deformation.

 Reprocessing, re-sterilisation, reuse and/or technical modification of the catheter by the user or operator!

Injury to blood vessels and intracardiac injuries, infection, sepsis, death of the patient

- The Vanguard Steerable Diagnostic Catheters SJA | SJF are disposable products and only intended for single use on one patient.
- The catheter may not be reprocessed, re-sterilised and reused by the user or operator.
- The catheter may not be technically modified.

 Devices emitting leakage currents are operated in the immediate vicinity of the patient!

Leakage currents on catheters, lethal arrhythmias

- Ensure that no devices with leakage currents are located and operated in the immediate vicinity of the patient.

-  The sheath used is too small for the catheter!
- The catheter is stuck, removal/replacement of the sheath and catheter/extension of surgery**
- The sheath used must be matched to the catheter before starting the procedure.
-  Pacing in the wrong place of the heart tissue!
- Ventricular arrhythmias requiring external defibrillation, possibly associated with skin burns**
- The medical device may only be used by a trained cardiologist specialised in electrophysiology.
-  Using the catheter for ablation or internal cardioversion!
- Perforation, arrhythmia, embolism, thrombus formation and/or death of the patient**
- Do not use the catheter for ablation or internal cardioversion, apply the catheter only in accordance with its intended use.
- The medical device may only be used by a trained cardiologist specialised in electrophysiology.
-  The catheter is stuck in the heart or vessel!
- Surgical intervention, repair of injured tissue and/or the damaged valve**
- Avoid twisting the catheter too much and ensure good imaging when positioning it around the chordae tendineae.
- The medical device may only be used by a trained cardiologist specialised in electrophysiology.
- Take particular care during use on patients with implanted cardiac electrodes.
-  Excessive mechanical stress on the catheter due to manual bending, twisting and curving!
- Impairment and/or failure of curvature control, patient injury**
- Avoid excessive mechanical stress on the catheter.
- The catheter must not be subjected to any mechanical stress and/or kinks by manual deformation.
-  Use of long introducer sheaths!
- Thromboembolic events**
- Pre-flush the lumen of the sheath with heparinised intravenous solution for infusion.
- Maintain the patency of the sheath's lumen.
-  Cleaning of the catheter with organic solvents e.g. alcohol before/during use on the patient!
- Toxic reactions, inflammation or allergic reactions**
- Do not treat/wipe the catheter with organic solvents.
- Only use water or sterile saline solution to wipe the catheter.
-  The electrodes of the Vanguard Steerable Diagnostic Catheter SJA | SJF are in contact with the ablation electrode!
- Damage to the catheter and/or charring, clots, embolism**
- When using the Vanguard Steerable Diagnostic Catheter SJA | SJF, ensure good imaging.
- All contact between the electrodes of the Vanguard Steerable Diagnostic Catheter SJA | SJF device and the electrodes of an ablation catheter needs to be avoided.
-  The EGM recorder is active when connected to the catheter!
- Triggering of arrhythmias**
- First connect the catheter to the EGM recorder. Only switch on the recorder once the catheter is connected.

 Use of connecting cables with unprotected plug contacts!

Death of the patient

- Only use cables with shielded plug contacts.
- The Vanguard Steerable Diagnostic Catheter SJA | SJF must be used in conjunction with compatible interconnect cables manufactured by *Abbott Medical (formerly Irvine Biomedical, Inc. a St. Jude Medical Company, 2375 Morse Avenue, Irvine, California 92614 USA)*.

DIRECTIONS FOR USE

It is the responsibility of the physician to use a suitable surgical procedure and technology. The procedure described in the instructions for use is for information purposes only. Each physician must apply, supplement or adapt the information in the instructions for use according to his/her medical training and clinical experience. The medical device may only be used by a trained cardiologist specialised in electrophysiology. Cardiac mapping procedures may only be performed by appropriately trained professionals in a regular and fully equipped electrophysiological (EP) laboratory.

The catheters are used in conjunction with compatible connection cables from *Abbott Medical (formerly Irvine Biomedical, Inc. a St. Jude Medical Company, 2375 Morse Avenue, Irvine, California 92614 USA)* and suitable devices. All Vanguard Steerable Diagnostic Catheters SJA | SJF are as compatible with accessories (e.g. cables) as the corresponding products of the original manufacturers *Abbott Medical (formerly St. Jude Medical, Inc.; refer to technical data)*.

Handling instructions

Handling the catheter requires great care to avoid injuries, perforations and tamponades to the heart and vascular system.

Storage

The product must be stored in a cool, dry and light-protected area. Storage temperatures range from 5 °C to 25 °C.

- Information about storage conditions is also provided on the label on the product packaging.

Safety measures before use

- Before using the catheter or external pacemaker, you must read and understand the respective instructions for use.
- Observe the relevant instructions for use when connecting and using external devices. Only use approved medical devices while observing the instructions for use.
- All static electricity must be kept away from the catheter system. You must also take particular care to ensure that the operating table and the electronic devices which are used (e.g. X-rays) are earthed sufficiently and centrally.
- It is a pre-requirement for intracardiac potential recording and/or stimulation that specialist medical personnel are trained sufficiently and that a fully equipped electrophysiology laboratory is present. The procedure itself may only be carried out by a trained cardiologist specialized in electrophysiology.

Removal from packaging

- Check that the packaging for integrity and the sterilization indicator (yellow = sterile). If the packaging is damaged or the sterilization indicator is not yellow, discard the product and use a new one for the procedure.
- Do not use the device once its "use by" date has expired.
- Do not use the device if its label is incomplete or illegible.
- Take the catheter out under aseptic conditions and transfer it to a sterile working environment.
- Ensure that the packaging material is disposed of safely.

- Check the catheter for any obvious defects such as kinks in the shaft, loose electrodes or damaged plugs and, if necessary, swap it for another catheter. The catheter must not be used if there is any obvious damage.

Application

- ✓ To avoid injury to the vascular system and the heart, the procedure should always be performed using appropriate visualisation techniques (e.g. fluoroscopy, ultrasound).
- ✓ An emergency power supply must be available to ensure continuous treatment.
- ✓ An external pacemaker, defibrillator, intubation kit, oxygen connection and emergency medication must be ready and available throughout the procedure.
- ✓ Defibrillation, emergency pacemaker therapy and acute cardioversion must be possible during the entire treatment.
- ✓ The principles of asepsis must be observed.
- ✓ In order to prevent delays in the process, a sterile substitute should be available during treatment.
- Under consideration of the aseptic rules, access a large central blood vessel through which you can push the catheter in.
- Use the appropriate connection cable to connect the catheter to a suitable EGM recorder/stimulation device.
- Switch on the EGM recorder after connecting it to the catheter.
- Before inserting and removing the catheter, make sure that its tip is in straight position.
- Advance the catheter using appropriate visualisation techniques and electrocardiogram and position it at the desired endocardial location.

NOTE: To facilitate positioning of the catheter tip on the Vanguard Steerable Diagnostic Catheter SJA| SJF models, it can be deflected using the thumb knob. To do this, the thumb knob is pushed forward (distal direction). If the thumb knob is pushed backwards, the tip straightens again. The catheter is straight when the thumb knob is pushed fully backwards (proximal end position).

NOTE: If resistance occurs, the catheter must not be pulled back with force. In this case, the catheter position must be checked by fluoroscopy examination.

NOTE: The catheters may only be used with approved recording devices and generators that are compatible with appropriate catheters made by *Abbott Medical (formerly St. Jude Medical, Inc.)*. Observe the information and Technical Manual of the equipment manufacturer. Use the cables supplied by *Abbott Medical (formerly St. Jude Medical, Inc.)* that are intended for use with the corresponding catheters made by the same company.

Safety measures after use

Following the clinical procedure, the catheter must be checked for intactness for patient safety purposes. Pay attention to breaks and fragmentation of the shaft and tip. All actively implanted devices of the patient must be checked for functionality.

Disposal

Any used products must be considered as contaminated hospital waste. The catheter and its packaging must be disposed of accordance with the applicable national laws and directives.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Current medical practice guidelines recommend catheter ablation as suitable therapeutic approach for a number of cardiac arrhythmias after weighing potential benefits and risks. Diagnostic catheters such as the Vanguard Steerable Diagnostic Catheters SJA | SJF are indispensable for conducting such ablation procedures.

Clinical data from literature and market surveillance demonstrate state of the art acute success rates for procedures performed using the Vanguard Steerable Diagnostic Catheters SJA | SJF in common application fields of catheter ablation.

Type and frequency of procedural complications likewise remain within the expected range according to current literature.

With that, the Vanguard Steerable Diagnostic Catheters SJA | SJF is a suitable and safe product for diagnosis of cardiac arrhythmias in the context of electrophysiological studies.

The Summary of Safety and Clinical Performance (SSCP), which includes clinical benefit, can be located as follows after the launch of European database on medical devices (Eudamed): <https://ec.europa.eu/tools/eudamed>. The Vanguard SJA | SJF catheters basic UDI-DI can be used to find the corresponding SSCP, BASIC UDI-DI 4038134EP-SJABCDEF00000NE

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©Vanguard AG 1x00003457_REV_C_15-11-2023