

CRITERIA FOR REVIEW OF INCOMING TRANSPORT - BOXES

The review of transport boxes will be under the following criteria, which can result to a rejection of the device

Criteria

- maximum reprocessing/full refurbishment amount has been reached
- kink in the shaft
- surface damaged
- time period until OEM Use Before Date is less than 13 months at goods receipt at Vanguard
- Label is missing
- product has been reprocessed/fully refurbished by third parties
- heavy decontamination/ staining
- product cannot be fully refurbished

The listed criteria may not complete.